

## Product datasheet for **HP212390**

### ARID1B Human qPCR Primer Pair (NM\_017519)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57492
Forward Sequence:	CCAGTCAACTGGCAGCAATTCC
Reverse Sequence:	CACTCACATCTGAGAATGGGTCG
ACCN:	<u>NM_017519</u> , <u>NM_017519.1</u> , <u>NM_017519.2</u> , <u>BC009543</u> , <u>BU733862</u>
UniProt ID:	<u>Q8NFD5</u>
Synonyms:	6A3-5; BAF250B; BRIGHT; CSS1; DAN15; ELD/OSA1; MRD12; OSA2; P250R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)