

## Product datasheet for **HP212414**

### INO80 Human qPCR Primer Pair (NM\_017553)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                             |
| Gene ID:          | 54617                                                                                                                                                                                                                                                                                                                                         |
| Forward Sequence: | CCACCAGAGATACCTGAGGAAC                                                                                                                                                                                                                                                                                                                        |
| Reverse Sequence: | CAACCTGGTCTGAGTAGCCACT                                                                                                                                                                                                                                                                                                                        |
| ACCN:             | <a href="#">BC012578</a> , <a href="#">NM_017553</a> , <a href="#">NM_017553.1</a> , <a href="#">NM_017553.2</a> , <a href="#">BC046115</a> , <a href="#">BC057299</a> , <a href="#">BC068054</a> , <a href="#">BC105993</a> , <a href="#">BC144681</a> , <a href="#">BC146785</a> , <a href="#">BX350023</a> , <a href="#">NM_017553.3</a>   |
| UniProt ID:       | <a href="#">Q9ULG1</a>                                                                                                                                                                                                                                                                                                                        |
| Synonyms:         | INO80A; INOC1                                                                                                                                                                                                                                                                                                                                 |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.                                                                                                                                                |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                               |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                  |



[View online »](#)