

## Product datasheet for **HP212441**

### BTG4 Human qPCR Primer Pair (NM\_017589)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54766
Forward Sequence:	GAAGTCACTGGCACTCTGATTGC
Reverse Sequence:	TGGTCATCTCCTTCGGAAGTCC
ACCN:	<a href="#">BC031045</a> , <a href="#">NM_017589</a> , <a href="#">NM_017589.2</a> , <a href="#">NM_017589.3</a> , <a href="#">BC031045.1</a> , <a href="#">BC022056</a> , <a href="#">BX099301</a> , <a href="#">NM_017589.4</a>
UniProt ID:	<a href="#">Q9NY30</a>
Synonyms:	APRO3; OOMD8; PC3B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)