

Product datasheet for **HP212450**

GOLGA2P5 Human qPCR Primer Pair (NM_017600)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55592
Forward Sequence:	CAGAAGGAACCTGAGGCAGCAG
Reverse Sequence:	TGGCTCTGGTACAGTCCGATGT
ACCN:	<u>NM_017600</u> , <u>NM_017600.1</u> , <u>BC034587</u> , <u>BC048272</u> , <u>BC092435</u> , <u>BC109389</u> , <u>BC110077</u> , <u>BC127181</u> , <u>BC141930</u> , <u>BQ366018</u> , <u>BX504274</u> , <u>BX647364</u>
Synonyms:	DKFZp434A0411; DKFZp434K1926; DKFZp434M0331; GOLGA2L1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)