

## Product datasheet for **HP212455**

### RNF111 Human qPCR Primer Pair (NM\_017610)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54778
Forward Sequence:	CTGTCACCAGTAGCCAACCTTC
Reverse Sequence:	GCCAGTTTCCGTTACAGCACTAG
ACCN:	<u><a href="#">NM_017610</a></u> , <u><a href="#">NM_017610.1</a></u> , <u><a href="#">NM_017610.2</a></u> , <u><a href="#">NM_017610.3</a></u> , <u><a href="#">NM_017610.4</a></u> , <u><a href="#">NM_017610.5</a></u> , <u><a href="#">NM_017610.6</a></u> , <u><a href="#">NM_017610.7</a></u> , <u><a href="#">BC060862</a></u> , <u><a href="#">BC010369</a></u> , <u><a href="#">BC020984</a></u> , <u><a href="#">BM986870</a></u> , <u><a href="#">BX538130</a></u> , <u><a href="#">BX647259</a></u> , <u><a href="#">NM_017610.8</a></u>
UniProt ID:	<u><a href="#">Q6ZNA4</a></u>
Synonyms:	ARK; hRNF111
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)