

## Product datasheet for **HP212471**

### DDX60 Human qPCR Primer Pair (NM\_017631)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55601
Forward Sequence:	GGTGTTTTACCAGGGAGTATCG
Reverse Sequence:	CCAGTTTTGGCGATGAGGAGCA
ACCN:	<a href="#">NM_017631</a> , <a href="#">NM_017631.1</a> , <a href="#">NM_017631.2</a> , <a href="#">NM_017631.3</a> , <a href="#">NM_017631.4</a> , <a href="#">NM_017631.5</a> , <a href="#">BC038115</a> , <a href="#">BC038115.2</a> , <a href="#">BC007820</a> , <a href="#">BC020601</a> , <a href="#">BC074781</a> , <a href="#">BX648655</a>
UniProt ID:	<a href="#">Q8IY21</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)