

Product datasheet for **HP212473**

TENT5A Human qPCR Primer Pair (NM_017633)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55603
Forward Sequence:	TGCCTGTTGGACTTCTTACCCG
Reverse Sequence:	GACAGGGATATAAGACTCCATCG
ACCN:	<u>NM_017633</u> , <u>NM_017633.1</u> , <u>NM_017633.2</u> , <u>BC000683</u> , <u>BC000683.2</u> , <u>BC007351</u> , <u>BX648876</u> , <u>NM_017633.3</u>
UniProt ID:	<u>Q96IP4</u>
Synonyms:	C6orf37; FAM46A; OI18; XTP11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)