

Product datasheet for **HP212476**

BNC2 Human qPCR Primer Pair (NM_017637)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54796
Forward Sequence:	ACTCTGCGGGACTATGTCCGAG
Reverse Sequence:	ACCGCAGAACTGCTGAAGGGT
ACCN:	<u>NM_017637</u> , <u>NM_017637.1</u> , <u>NM_017637.2</u> , <u>NM_017637.3</u> , <u>NM_017637.4</u> , <u>NM_017637.5</u> , <u>BC037160</u> , <u>BC156477</u> , <u>BC172513</u> , <u>BX641030</u> , <u>NM_017637.6</u>
UniProt ID:	<u>Q6ZN30</u>
Synonyms:	BSN2; LUTO
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)