

Product datasheet for **HP212488**

AHI1 Human qPCR Primer Pair (NM_017651)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54806
Forward Sequence:	GCTCAGTAGACACAGAACCTGG
Reverse Sequence:	CTCCTGCATTTAGTGAGAAGAGG
ACCN:	NM_017651 , NM_017651.1 , NM_017651.2 , NM_017651.3 , NM_017651.4 , BC015019 , BC029417 , BC039559 , BC050607 , BC065712 , BC094800 , BX643029
UniProt ID:	Q8N157
Synonyms:	AHI-1; dj71N10.1; JBTS3; ORF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)