

Product datasheet for **HP212490**

SAMD9 Human qPCR Primer Pair (NM_017654)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54809
Forward Sequence:	GGGAACCTACCTTGGCTATGCAC
Reverse Sequence:	CGTATTCCTGACGGTTCATTGCC
ACCN:	<u>NM_017654</u> , <u>NM_017654.1</u> , <u>NM_017654.2</u> , <u>NM_017654.3</u> , <u>BC065769</u> , <u>BC132773</u> , <u>BC132775</u> , <u>BC150249</u> , <u>BF109390</u> , <u>BG615970</u> , <u>BX647072</u> , <u>NM_017654.4</u>
UniProt ID:	<u>Q5K651</u>
Synonyms:	C7orf5; DRIF1; M7MLS2; MIRAGE; NFTC; OEF1; OEF2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)