

Product datasheet for **HP212500**

VPS50 Human qPCR Primer Pair (NM_017667)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55610
Forward Sequence:	TCCTTCAGCCACATCTGGATGC
Reverse Sequence:	AGTTCAGTGGCGTTGAGACTG
ACCN:	BC017888 , NM_017667 , NM_017667.1 , NM_017667.2 , NM_017667.3 , BC035984 , BC047757 , BC108708 , BC132740 , BC132742 , BC152424 , BC160012 , BE671517 , NM_017667.4
UniProt ID:	Q96JG6
Synonyms:	CCDC132; VPS54L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)