

Product datasheet for **HP212530**

ASAP3 Human qPCR Primer Pair (NM_017707)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55616
Forward Sequence:	CGGAAGAACTCAGGATGTGGCT
Reverse Sequence:	AGGCAGCCATACTTGACTCCAC
ACCN:	<u>NM_017707</u> , <u>NM_017707.1</u> , <u>NM_017707.2</u> , <u>NM_017707.3</u> , <u>BC060786</u> , <u>BC023519</u> , <u>BC052305</u> , <u>NM_017707.4</u>
UniProt ID:	<u>Q8TDY4</u>
Synonyms:	ACAP4; CENTB6; DDEFL1; UPLC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)