

## Product datasheet for **HP212555**

### ZCCHC2 Human qPCR Primer Pair (NM\_017742)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54877
Forward Sequence:	AGACTCTCCACTCCTCAGCATG
Reverse Sequence:	GGACTAGATGGAGAACTGTAACTC
ACCN:	<a href="#">BC006340</a> , <a href="#">NM_017742</a> , <a href="#">NM_017742.1</a> , <a href="#">NM_017742.2</a> , <a href="#">NM_017742.4</a> , <a href="#">NM_017742.5</a> , <a href="#">BC006340.1</a> , <a href="#">BC047606</a> , <a href="#">BC137435</a> , <a href="#">BI023990</a> , <a href="#">BQ001246</a> , <a href="#">BU675299</a> , <a href="#">BX485381</a>
UniProt ID:	<a href="#">Q9C0B9</a>
Synonyms:	C18orf49
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)