

Product datasheet for **HP212567**

NSUN2 Human qPCR Primer Pair (NM_017755)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54888
Forward Sequence:	ACCTGGCTCAAAGACCACACAG
Reverse Sequence:	TGGCTTGATGGACGAGCAGGTA
ACCN:	<u>BC001041</u> , <u>NM_017755</u> , <u>NM_017755.1</u> , <u>NM_017755.2</u> , <u>NM_017755.3</u> , <u>NM_017755.4</u> , <u>NM_017755.5</u> , <u>BC137083</u> , <u>BF000644</u> , <u>BG387061</u> , <u>BQ689094</u> , <u>NM_017755.6</u>
UniProt ID:	<u>Q08J23</u>
Synonyms:	MISU; MRT5; SAKI; TRM4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)