

Product datasheet for **HP212584**

TTC19 Human qPCR Primer Pair (NM_017775)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54902
Forward Sequence:	CGCTACCTTCTGTTCTCCAAGC
Reverse Sequence:	AGTAGCCAGGTCACTCATCAGC
ACCN:	<u>NM_017775</u> , <u>NM_017775.1</u> , <u>NM_017775.2</u> , <u>NM_017775.3</u> , <u>BC105128</u> , <u>BC105128.1</u> , <u>BC003058</u> , <u>BC011698</u> , <u>BC073796</u> , <u>BC112107</u>
UniProt ID:	<u>Q6DKK2</u>
Synonyms:	2010204O13Rik; MC3DN2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)