

Product datasheet for **HP212599**

BTBD2 Human qPCR Primer Pair (NM_017797)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55643
Forward Sequence:	AAGGAGTGCAGCATCAACCGCT
Reverse Sequence:	CACCACGAAGATGCGCTTGTG
ACCN:	<u>NM_017797</u> , <u>NM_017797.1</u> , <u>NM_017797.2</u> , <u>NM_017797.3</u> , <u>BC000564</u> , <u>BC008035</u> , <u>BC053986</u> , <u>BC069199</u> , <u>BC094820</u> , <u>BC152718</u> , <u>BC156974</u> , <u>BU849231</u> , <u>NM_017797.4</u>
UniProt ID:	<u>Q9BX70</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)