

Product datasheet for **HP212603**

HEATR2 (DNAAF5) Human qPCR Primer Pair (NM_017802)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54919
Forward Sequence:	CTCTTGGACGTGCTGCTGACAA
Reverse Sequence:	GGACCAATGTGCTTGCGGTAGA
ACCN:	NM_017802 , NM_017802.1 , NM_017802.2 , NM_017802.3 , BC047240 , BC047240.1 , BC010850 , BC072425 , BC094840 , BM743055 , NM_017802.4
UniProt ID:	Q86Y56
Synonyms:	CILD18; HEATR2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)