

Product datasheet for **HP212605**

RAIN (RASIP1) Human qPCR Primer Pair (NM_017805)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54922
Forward Sequence:	CAGGAGAAGGTGCTGGAATGG
Reverse Sequence:	GACAGACTGCTGGAAGGTACAC
ACCN:	<u>NM_017805</u> , <u>NM_017805.1</u> , <u>NM_017805.2</u> , <u>BC028614</u> , <u>BC028614.1</u> , <u>BC021860</u> , <u>BC042111</u> , <u>BX389914</u>
UniProt ID:	<u>Q5U651</u>
Synonyms:	RAIN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)