

Product datasheet for **HP212614**

HAUS4 Human qPCR Primer Pair (NM_017815)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54930
Forward Sequence:	GCAGATGTTGCTGCTGGAGAAG
Reverse Sequence:	ACCGCACTTGACTCCAGGTAC
ACCN:	NM_017815 , NM_017815.1 , NM_017815.2 , BC002554 , BC002554.2 , BC001916 , BQ066232 , BX247981 , BX247995 , BX248256 , NM_017815.3
UniProt ID:	Q9H6D7
Synonyms:	C14orf94
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)