

Product datasheet for **HP212630**

C9orf6 (FAM206A) Human qPCR Primer Pair (NM_017832)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54942
Forward Sequence:	GATCAGTACCAACTGTAGCAGAC
Reverse Sequence:	TCTTCACCATCTGAGCAGTAAATC
ACCN:	<u>NM_017832</u> , <u>NM_017832.1</u> , <u>NM_017832.2</u> , <u>NM_017832.3</u> , <u>BC015795</u> , <u>BC015795.1</u> , <u>BC051827</u> , <u>NM_017832.4</u>
UniProt ID:	<u>Q9NX38</u>
Synonyms:	C9orf6; CG-8; FAM206A; Simiate
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)