

Product datasheet for **HP212639**

TRSPAP1 (TRNAU1AP) Human qPCR Primer Pair (NM_017846)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54952
Forward Sequence:	TGGCTATACCCAGAGCACCATG
Reverse Sequence:	AGTGACAGTCCATCAGAGCGTC
ACCN:	NM_017846 , NM_017846.1 , NM_017846.2 , NM_017846.3 , NM_017846.4 , BC000680 , BC000680.2 , BC035855 , BC039879 , BI602122 , BQ014404 , NM_017846.5
UniProt ID:	Q9NX07
Synonyms:	PRO1902; SECP43; TRSPAP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)