

Product datasheet for **HP212644**

PARP16 Human qPCR Primer Pair (NM_017851)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 54956 |
| Forward Sequence: | GCATTTTCATGGTAGCCGCCTAG |
| Reverse Sequence: | CAAGTCACTGGTGAGGTAGGTC |
| ACCN: | NM_017851 , NM_017851.1 , NM_017851.2 , NM_017851.3 , NM_017851.4 , NM_017851.5 , BC006389 , BC006389.1 , BC031074 , NM_017851.6 |
| UniProt ID: | Q8N5Y8 |
| Synonyms: | ARTD15; C15orf30; pART15 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)