

Product datasheet for **HP212649**

SSH3 Human qPCR Primer Pair (NM_017857)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54961
Forward Sequence:	GCCTCAATGAGTGGACGGCTAT
Reverse Sequence:	ACTCTCCAGGTCCTGACATCC
ACCN:	<u>NM_017857</u> , <u>NM_017857.1</u> , <u>NM_017857.2</u> , <u>NM_017857.3</u> , <u>BC007709</u> , <u>BC007709.1</u> , <u>BC004176</u> , <u>BC004210</u> , <u>BU689839</u> , <u>NM_017857.4</u>
UniProt ID:	<u>Q8TE77</u>
Synonyms:	SSH3L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)