

Product datasheet for **HP212660**

INT11 (CPSF3L) Human qPCR Primer Pair (NM_017871)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54973
Forward Sequence:	GATCCGCAAGACTTTCGTGCAG
Reverse Sequence:	AGCATTCTGCGTGGCAAACA
ACCN:	<u>BC008041</u> , <u>NM_017871</u> , <u>NM_017871.1</u> , <u>NM_017871.2</u> , <u>NM_017871.3</u> , <u>NM_017871.4</u> , <u>NM_017871.5</u> , <u>BC007978</u> , <u>BC007978.2</u> , <u>BC000675</u> , <u>BC013904</u> , <u>BC020199</u> , <u>BC042444</u> , <u>BC047576</u> , <u>BK005673</u> , <u>BK005728</u>
UniProt ID:	<u>Q5TA45</u>
Synonyms:	CPSF3L; CPSF73L; INT11; RC-68; RC68
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)