

## Product datasheet for **HP212695**

### **C12orf48 (PARPBP) Human qPCR Primer Pair (NM\_017915)**

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                             |
| Gene ID:          | 55010                                                                                                                                                                                                                                                                                                                                         |
| Forward Sequence: | CTTAGCACCCTGACATCAGTCC                                                                                                                                                                                                                                                                                                                        |
| Reverse Sequence: | GGAGCATTAGCTGCTTCTTCGTC                                                                                                                                                                                                                                                                                                                       |
| ACCN:             | <u>NM_017915</u> , <u>NM_017915.1</u> , <u>NM_017915.2</u> , <u>NM_017915.3</u> , <u>NM_017915.4</u> , <u>BC014313</u> , <u>BC018903</u> , <u>BC030962</u> , <u>BC050696</u> , <u>BC070270</u> , <u>BC098297</u> , <u>BC098313</u> , <u>BC099734</u> , <u>BC099844</u> , <u>BG285036</u> , <u>NM_017915.5</u>                                 |
| UniProt ID:       | <u>Q9NWS1</u>                                                                                                                                                                                                                                                                                                                                 |
| Synonyms:         | AROM; C12orf48; PARI                                                                                                                                                                                                                                                                                                                          |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.                                                                                                                                                |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                               |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                  |



[View online »](#)