

Product datasheet for **HP212717**

USP47 Human qPCR Primer Pair (NM_017944)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55031
Forward Sequence:	CTGCTGTTGGAGACGAGAAAGC
Reverse Sequence:	GCACGAACAGTTATAGGAGCAGC
ACCN:	<u>NM_017944</u> , <u>NM_017944.1</u> , <u>NM_017944.2</u> , <u>NM_017944.3</u> , <u>BC000226</u> , <u>BC017795</u> , <u>BC030562</u> , <u>BC047044</u> , <u>BC071559</u> , <u>BC148587</u> , <u>BG752249</u> , <u>BM314303</u> , <u>NM_017944.4</u>
UniProt ID:	<u>Q96K76</u>
Synonyms:	TRFP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)