

Product datasheet for **HP212727**

TRMT12 Human qPCR Primer Pair (NM_017956)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55039
Forward Sequence:	CTCGAACTCCAGCAGTGACACT
Reverse Sequence:	AGCACAGGACAACGATGCCACT
ACCN:	<u>NM_017956</u> , <u>NM_017956.1</u> , <u>NM_017956.2</u> , <u>NM_017956.3</u> , <u>BC011713</u> , <u>BC011713.2</u> , <u>BC003057</u> , <u>BP342608</u> , <u>NM_017956.4</u>
UniProt ID:	<u>Q53H54</u>
Synonyms:	TRM12; TYW2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)