

Product datasheet for **HP212830**

NETO2 Human qPCR Primer Pair (NM_018092)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81831
Forward Sequence:	GTAGTCGGCTTAGCAGGTTTCG
Reverse Sequence:	GGTATGCACAATTTGGACACCAT
ACCN:	NM_018092 , NM_018092.1 , NM_018092.2 , NM_018092.3 , NM_018092.4 , BC012381 , BC065504 , BC098380 , BC108727 , BC152765 , BC156955 , BX538063 , NM_018092.5
UniProt ID:	Q8NC67
Synonyms:	BTCL2; NEOT2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)