

Product datasheet for **HP212845**

TMEM38B Human qPCR Primer Pair (NM_018112)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55151
Forward Sequence:	GTCAGTGATGGCGGTGAAACGT
Reverse Sequence:	TCAATGGAGGCTCTGCAAGCAG
ACCN:	<u>NM_018112</u> , <u>NM_018112.1</u> , <u>NM_018112.2</u> , <u>BC000049</u> , <u>BC031938</u> , <u>BU728078</u> , <u>BX112234</u> , <u>NM_018112.3</u>
UniProt ID:	<u>Q9NVV0</u>
Synonyms:	bA219P18.1; C9orf87; D4Ert89e; OI14; TRIC-B; TRICB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)