

Product datasheet for **HP212860**

SHQ1 Human qPCR Primer Pair (NM_018130)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55164
Forward Sequence:	CAGGACAGCAAACCACTTGTGG
Reverse Sequence:	TCCATAGTGCCACTGCGGATTC
ACCN:	NM_018130 , NM_018130.1 , NM_018130.2 , BC032671 , BC032671.1 , BC017204 , BC017274 , BC025270 , BC039830 , BC047879 , BX641074 , NM_018130.3
UniProt ID:	Q6PI26
Synonyms:	GRIM-1; Shq1p
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)