

Product datasheet for **HP212914**

MACO1 Human qPCR Primer Pair (NM_018202)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55219
Forward Sequence:	CTGAGCAGGAAGCCCGAAGTTT
Reverse Sequence:	CGATTCCGTAAGGTTTCGGTGC
ACCN:	NM_018202 , NM_018202.1 , NM_018202.2 , NM_018202.3 , NM_018202.4 , NM_018202.5 , BC032427 , BC032427.1 , BC007687 , BC020801 , NM_018202.6
UniProt ID:	Q8N5G2
Synonyms:	MACOILIN; TMEM57
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)