

Product datasheet for **HP212928**

MOBK1B (MOB1A) Human qPCR Primer Pair (NM_018221)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55233
Forward Sequence:	TGTTGCCTGAGGGAGAGGATCT
Reverse Sequence:	GCAGACATGACTGGACAGCTTG
ACCN:	<u>NM_018221</u> , <u>NM_018221.1</u> , <u>NM_018221.2</u> , <u>NM_018221.3</u> , <u>NM_018221.4</u> , <u>BC003398</u> , <u>BC003398.1</u> , <u>BI825399</u> , <u>BU172616</u> , <u>NM_018221.5</u>
UniProt ID:	<u>Q9H8S9</u>
Synonyms:	C2orf6; MATS1; MOB1; Mob4B; MOBK1B; MOBKL1B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)