

Product datasheet for **HP212949**

Transmembrane protein 30A (TMEM30A) Human qPCR Primer Pair (NM_018247)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55754
Forward Sequence:	TGCTCCTGTGGAGCTATTGCC
Reverse Sequence:	TCTGTCCACCAAGCAATACCTTTC
ACCN:	<u>NM_018247</u> , <u>NM_018247.1</u> , <u>NM_018247.2</u> , <u>NM_018247.3</u> , <u>BC009006</u> , <u>NM_018247.4</u>
UniProt ID:	<u>Q9NV96</u>
Synonyms:	C6orf67; CDC50A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)