

Product datasheet for **HP212950**

NEI3 (NEIL3) Human qPCR Primer Pair (NM_018248)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55247
Forward Sequence:	AGTGGTCTCCACCCAGCTGTTA
Reverse Sequence:	AGAGCAAGTCCTGCTTTACGGC
ACCN:	<u>NM_018248</u> , <u>NM_018248.1</u> , <u>NM_018248.2</u> , <u>BC025954</u> , <u>BC025954.1</u> , <u>BE865481</u> , <u>BP380910</u> , <u>NM_018248.3</u>
UniProt ID:	<u>Q8TAT5</u>
Synonyms:	FPG2; FPG2; hFPG2; hNEI3; NEI3; ZGRF3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)