

Product datasheet for **HP212965**

WDR41 Human qPCR Primer Pair (NM_018268)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55255
Forward Sequence:	GTTGCCAGACTTCCAACAGGC
Reverse Sequence:	ACCTGTTGGTACAGGCTCAGCT
ACCN:	NM_018268 , NM_018268.2 , NM_018268.3 , BC040241 , BC040241.1 , BC001807 , BG113002 , BM818066 , BP276774 , BQ000743 , NM_018268.4
UniProt ID:	Q9HAD4
Synonyms:	MSTP048
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)