

Product datasheet for **HP212967**

C20orf20 (MRGBP) Human qPCR Primer Pair (NM_018270)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55257
Forward Sequence:	TCGGGACAAGTTCAGCCAGAAC
Reverse Sequence:	TCCTCTCTGGATTCGGGAATGG
ACCN:	<u>NM_018270</u> , <u>NM_018270.1</u> , <u>NM_018270.2</u> , <u>NM_018270.3</u> , <u>NM_018270.4</u> , <u>NM_018270.5</u> , <u>BC018841</u> , <u>BC018841.2</u> , <u>BC009889</u> , <u>BC014235</u> , <u>BM920599</u> , <u>BQ445915</u> , <u>BU902487</u> , <u>NM_018270.6</u>
UniProt ID:	<u>Q9NV56</u>
Synonyms:	C20orf20; Eaf7; MRG15BP; URCC4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)