

Product datasheet for **HP212970**

TRAPPC14 Human qPCR Primer Pair (NM_018275)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55262
Forward Sequence:	CCTTCACTCAGAGCATCTACACC
Reverse Sequence:	GTGAAACGCTCGTAGGTCCGAA
ACCN:	NM_018275 , NM_018275.2 , NM_018275.3 , NM_018275.4 , BC015722 , BC015722.2 , BC002749 , BC094727 , NM_018275.5
UniProt ID:	Q8WVR3
Synonyms:	C7orf43; MCPH25; TRAPPC14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)