

Product datasheet for **HP212975**

ECHDC2 Human qPCR Primer Pair (NM_018281)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55268
Forward Sequence:	TGGAGTGAAGGGCGTGTCTGT
Reverse Sequence:	CGTGGTCTCAATCAGTCCCATG
ACCN:	NM_018281 , NM_018281.1 , NM_018281.2 , NM_018281.3 , BC044574 , BC038389 , BC040480 , BC051887 , BX647186 , BX648522 , NM_018281.4
UniProt ID:	Q86YB7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)