

## Product datasheet for **HP212977**

### GBP3 Human qPCR Primer Pair (NM\_018284)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2635
Forward Sequence:	TCCAAGGAGTCTGTGACCGATG
Reverse Sequence:	CCTCCACCATTTTGCTGAAGCC
ACCN:	<a href="#">NM_018284</a> , <a href="#">NM_018284.1</a> , <a href="#">NM_018284.2</a> , <a href="#">BC017992</a> , <a href="#">BC063819</a> , <a href="#">BC140837</a> , <a href="#">BC144650</a> , <a href="#">BP281050</a> , <a href="#">BU678920</a> , <a href="#">NM_018284.3</a>
UniProt ID:	<a href="#">Q9H0R5</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)