

Product datasheet for **HP212988**

NGLY1 Human qPCR Primer Pair (NM_018297)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55768
Forward Sequence:	CAGCACAAGTAGAACTGACAGGC
Reverse Sequence:	GTGTTGCCAAGCGACATCACCA
ACCN:	NM_018297 , NM_018297.1 , NM_018297.2 , NM_018297.3 , BC007226 , BC007226.1 , BC017220 , BC000963 , BQ017148 , NM_018297.4
UniProt ID:	Q96IV0
Synonyms:	CDDG; CDG1V; PNG-1; PNG1; PNGase
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)