

Product datasheet for **HP212994**

TBC1D23 Human qPCR Primer Pair (NM_018309)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55773
Forward Sequence:	CTTTGTGGTGGATTGCCGTCCT
Reverse Sequence:	GCCAGACTCAATGGACTGCTTC
ACCN:	NM_018309 , NM_018309.1 , NM_018309.2 , NM_018309.3 , NM_018309.4 , BC010221 , BC020955 , BX648032 , NM_018309.5
UniProt ID:	Q9NUY8
Synonyms:	NS4ATP1; PCH11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)