

Product datasheet for **HP213018**

SLC29A3 Human qPCR Primer Pair (NM_018344)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55315
Forward Sequence:	GGCAGTCTACTGCCATGGAAC
Reverse Sequence:	CAACGGCAAGGTAGCTCTCAAAG
ACCN:	NM_018344 , NM_018344.1 , NM_018344.2 , NM_018344.3 , NM_018344.4 , NM_018344.5 , BC120996 , BC000223 , BC041575 , BC050589 , BC063019 , BC120997
UniProt ID:	Q9BZD2
Synonyms:	ENT3; HCLAP; HJCD; PHID
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)