

Product datasheet for **HP213038**

Phytoceramidase (ACER3) Human qPCR Primer Pair (NM_018367)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55331
Forward Sequence:	CCTTGTTCTATCTTCACATCC
Reverse Sequence:	GAGGCTCAACAGGATCACTGG
ACCN:	NM_018367 , NM_018367.1 , NM_018367.2 , NM_018367.3 , NM_018367.4 , NM_018367.5 , NM_018367.6 , BC073853 , BC073853.1 , BC049837 , BC063034 , BG702017
UniProt ID:	Q9NUN7
Synonyms:	APHC; PHCA; PLDECO
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)