

Product datasheet for **HP213058**

C4orf21 (ZGRF1) Human qPCR Primer Pair (NM_018392)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55345
Forward Sequence:	TGGATGTCAGCCCTCTGGCTTA
Reverse Sequence:	CTTAGCCTCATGTGATGCAGCTG
ACCN:	<u>BC028584</u> , <u>NM_018392</u> , <u>NM_018392.1</u> , <u>NM_018392.2</u> , <u>NM_018392.3</u> , <u>NM_018392.4</u> , <u>BC044799</u> , <u>BC044799.1</u> , <u>BF056790</u> , <u>BM721145</u> , <u>BX647594</u> , <u>NM_018392.5</u>
UniProt ID:	<u>Q86YA3</u>
Synonyms:	C4orf21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)