

Product datasheet for **HP213069**

LAPTM4B Human qPCR Primer Pair (NM_018407)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55353
Forward Sequence:	GCTGTGTTTGGAACTGCTACCG
Reverse Sequence:	GCAGCACCATTACAGTGGCAT
ACCN:	<u>NM_018407</u> , <u>NM_018407.1</u> , <u>NM_018407.2</u> , <u>NM_018407.3</u> , <u>NM_018407.4</u> , <u>BC038117</u> , <u>BC038117.1</u> , <u>BC014129</u> , <u>BC031021</u>
UniProt ID:	<u>Q86VI4</u>
Synonyms:	LAPTM4beta; LC27
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)