

Product datasheet for **HP213086**

RRN3 Human qPCR Primer Pair (NM_018427)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54700
Forward Sequence:	CGGCAGGGTATTGAAGATGCTG
Reverse Sequence:	ACAGGATGCACCATCTGGTCGA
ACCN:	<u>NM_018427</u> , <u>NM_018427.1</u> , <u>NM_018427.2</u> , <u>NM_018427.3</u> , <u>NM_018427.4</u> , <u>BC132688</u> , <u>BC036182</u> , <u>BC071868</u> , <u>BC104660</u> , <u>BC132690</u> , <u>BU676519</u> , <u>NM_018427.5</u>
UniProt ID:	<u>Q9NYV6</u>
Synonyms:	A-270G1.2; TIFIA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)