

Product datasheet for **HP213090**

DOK5 Human qPCR Primer Pair (NM_018431)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55816
Forward Sequence:	TGACGATACCTCCAAGACTTTTGC
Reverse Sequence:	GCCAGTAAGTCAGGCTCTCCAA
ACCN:	<u>NM_018431</u> , <u>NM_018431.1</u> , <u>NM_018431.2</u> , <u>NM_018431.3</u> , <u>NM_018431.4</u> , <u>BC008992</u> , <u>BG212868</u> , <u>BI758267</u> , <u>BI916998</u> , <u>NM_018431.5</u>
UniProt ID:	<u>Q9P104</u>
Synonyms:	C20orf180; IRS-6; IRS6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)