

## Product datasheet for **HP213099**

### PDP1 Human qPCR Primer Pair (NM\_018444)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54704
Forward Sequence:	TTCTGGAGCCACTGCTTGTGTG
Reverse Sequence:	ACAGCGTGACTGCTGACCATGA
ACCN:	<a href="#">NM_018444</a> , <a href="#">NM_018444.1</a> , <a href="#">NM_018444.2</a> , <a href="#">NM_018444.3</a> , <a href="#">BC047619</a> , <a href="#">BC047619.1</a> , <a href="#">BC064978</a> , <a href="#">BC098343</a> , <a href="#">BM973633</a> , <a href="#">BQ422289</a> , <a href="#">BX099681</a> , <a href="#">NM_018444.4</a>
UniProt ID:	<a href="#">Q9P0J1</a>
Synonyms:	PDH; PDP; PDPC; PPM2A; PPM2C
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)