

Product datasheet for **HP213102**

TIP120A (CAND1) Human qPCR Primer Pair (NM_018448)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55832
Forward Sequence:	GCACTGTGAGTGTGCAGAGGAA
Reverse Sequence:	TTCACAGCCGTAACCACTGAGC
ACCN:	NM_018448 , NM_018448.1 , NM_018448.2 , NM_018448.3 , NM_018448.4 , BC050341 , BC050341.1 , BC004232 , BC026220 , BQ024785 , NM_018448.5
UniProt ID:	Q86VP6
Synonyms:	TIP120; TIP120A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)