

Product datasheet for **HP213112**

CISD1 Human qPCR Primer Pair (NM_018464)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55847
Forward Sequence:	CCTTCACATCCAGAAAGACAACC
Reverse Sequence:	CTCTTCGTTATGTTTTGTGTGAGC
ACCN:	NM_018464 , NM_018464.1 , NM_018464.2 , NM_018464.3 , NM_018464.4 , BC059168 , BC059168.1 , BC005962 , BC007043 , BC008474 , NM_018464.5
UniProt ID:	Q9NZ45
Synonyms:	C10orf70; MDS029; mitoNEET; ZCD1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)